



Bedinger Lab *In-vitro* liquid tomato pollen germination media protocol

Stock recipe

PGM	Stock Concentration	To make 5ml	10ml	20ml	30ml	40ml
PEG 4000	40%	3	6	12	18	24
Boric Acid	0.1%	0.5	1	2	3	4
Sucrose	40%	0.25	0.5	1	1.5	2
HEPES buffer	0.5M, pH6.0	0.2	0.4	0.8	1.2	1.6
Ca(NO ₃) ₂ ·4H ₂ O	0.1M	0.15	0.3	0.6	0.9	1.2
MgSO ₄ ·7H ₂ O	2%	0.05	0.1	0.2	0.3	0.4
KNO ₃	1% or 0.1M	0.05	0.1	0.2	0.3	0.4
H ₂ O	Double Distilled	0.8	1.6	3.2	4.8	5.6

1. Collect pollen from open and dehiscent flowers off a young tomato plant into 1.5 ml tubes with 1ml pollen germination media (PGM). Label the tube.
2. Mix well by inversion and minivortexing with a pollen buzzer. Record time as t_0 .
3. Dilute pollen/PGM solution 1:10 into a total of 300 μ l of new PGM.
 - a. I take 30 μ l of the pollen/PGM solution and add it to 270 μ l of PGM.
 - b. Remix by vortexing.
 - c. Use a cut tip to transfer.
4. After mixing aliquot 100 μ l of the diluted pollen/PGM into an individual well of a 96-well plate.
 - a. We view on an inverted scope because the pollen settles to the bottom of the well and is easily viewed. For an upright scope you'll have to use slides or a shallow well plate.

Notes:

- The quality of the PEG can be an issue. If available try batches from a couple different powder stocks. Getting this into solution does take time.
- All my stock solutions are filter sterilized; PEG and Sucrose take the most time to filter because at 40% they are viscous.
- All stock solutions are kept at -20 C and the media is made fresh each time. It wouldn't be a problem to make some media and keep it stored frozen as well and use it when you need to.